


Product Datasheet

Human ANGPTL3 ELISA Kit (orb215874)

Description	Human ANGPTL3 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S17-E460
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is around 33 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	53637 MW
Uniprot ID	Q9Y5C1
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
ANGPTL3
PicoKine
ELISA Kit
standar...