

## Product Datasheet

### CXCL5 ELISA Kit (Rat) (orb2159207)

<b>Description</b>	CXCL5 ELISA Kit (Rat)
<b>Range</b>	15.6 pg/ml - 1,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: A38-Q130
<b>Storage</b>	Store at 4C for 6 months, at -2°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernates, serum and plasma (heparin, EDTA, citrate)
<b>Uniprot ID</b>	<b>P97885</b>
<b>Sensitivity</b>	1 pg/ml
<b>NCBI</b>	<b>NP_071550.1</b>
<b>Expiration Date</b>	Please enquire.