


---

## Product Datasheet

### FGFR4 Antibody (orb215947)

Description	FGFR4 Antibody	100KD -
<b>Species/Host</b>	Rabbit	55KD -
<b>Reactivity</b>	Human	35KD -
<b>Conjugation</b>	Unconjugated	25KD -
<b>Tested Applications</b>	IHC, WB	15KD -
<b>Immunogen</b>	E.coli-derived human FGFR4 recombinant protein (Position: L22-H206). Human FGFR4 shares 86% and 84% amino acid (aa) sequences identity with mouse and rat FGFR4, respectively.	WB analysis of FGFR4 using anti-FGFR4 an...
<b>Form/Appearance</b>	Lyophilized	130KD - 100KD - 70KD - 55KD - 35KD - 25KD -
<b>Concentration</b>	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	WB analysis of FGFR4 using anti-FGFR4 an...
<b>Storage</b>	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
<b>Note</b>	For research use only	IHC analysis of FGFR4 using anti-FGFR4 a...
<b>Application notes</b>	WB: The detection limit for FGFR4 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
<b>Isotype</b>	Rabbit IgG	
<b>Clonality</b>	Polyclonal	
<b>MW</b>	87954 MW	
<b>Uniprot ID</b>	<b>P22455</b>	
<b>Expiration Date</b>	12 months from date of receipt.	