

Product Datasheet

CD40L ELISA Kit (Human) (orb2159506)

Description	CD40L ELISA Kit (Human)
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E108-L261
Storage	Store at 4C for 6 months, at -2°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum and plasma (heparin, EDTA)
Uniprot ID	P29965
Sensitivity	15 pg/ml
NCBI	NP_000065.1
Expiration Date	Please enquire.