

Product Datasheet

IGFBP-3 ELISA Kit (Human) [Out of stock] (orb2159579)

Description	IGFBP-3 ELISA Kit (Human) [Out of stock]
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A28Q-291
Storage	Store at 4C for 6 months, at -2°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum and plasma (heparin, EDTA)
Uniprot ID	P17936
Sensitivity	10 pg/ml
NCBI	NP_000589.2
Expiration Date	Please enquire.