

Product Datasheet

EPHB1/2/3/4 Colorimetric Cell-Based ELISA Kit (orb2160349)

Description	EPHB1/2/3/4 Colorimetric Cell-Based ELISA Kit
Range	At least 5000 cells/well
Tested Applications	ELISA
Preservatives	96-Well Microplate
Form/Appearance	96-Well Microplate
Storage	2°C to 8°C
Note	For research use only
Application notes	Assay Info: Assay Type: Cell-Based Subtype: None Detection Method: Colorimetric 450 nm
Uniprot ID	P54762 , P54760 , P29323 , P54753
Expiration Date	Please enquire.