

Product Datasheet

Human IgG Fab antibody (orb216255)

Description	Mouse monoclonal antibody to Human IgG Fab
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IH, WB
Immunogen	Human IgG Fab
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.05% Sodium Azide
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	ELISA: 1:5000-1:20000, WB: 1:5000-1:20000, IHC-P: 1:500-1:1000
Expiration Date	12 months from date of receipt.