

Product Datasheet

Guinea Pig IgG (H&L) antibody (HRP) (orb216287)

Description	Rabbit polyclonal antibody against Guinea Pig IgG (H&L) conjugated to HRP
Species/Host	Rabbit
Reactivity	Guinea pig
Conjugation	HRP
Tested Applications	ELISA, IH, WB
Immunogen	Guinea Pig IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Proclin
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	ELISA: 1:5000-1:20000, WB: 1:5000-1:20000, IHC-P: 1:500-1:1000
Expiration Date	12 months from date of receipt.