

## Product Datasheet

# Guinea Pig IgG (H&L) antibody (Biotin) (orb216327)

<b>Description</b>	Rabbit polyclonal antibody against Guinea Pig IgG (H&L) conjugated to Biotin
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Guinea pig
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IH, WB
<b>Immunogen</b>	Guinea Pig IgG
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Dilution Range</b>	ELISA: 1:5000-1:20000, WB: 1:5000-1:20000, IHC-P: 1:500-1:1000
<b>Expiration Date</b>	12 months from date of receipt.