

Product Datasheet

Guinea Pig IgG (H&L) antibody (FITC) (orb216367)

Description	Rabbit polyclonal antibody against Guinea Pig IgG (H&L) conjugated to FITC
Species/Host	Rabbit
Reactivity	Guinea pig
Conjugation	FITC
Tested Applications	ELISA, FC, IF, IH
Immunogen	Guinea Pig IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	IF/ICC: 1:250-1:2000, FACS: 1:250-1:2000
Expiration Date	12 months from date of receipt.