

Product Datasheet

Bradykinin Potentiator C peptide (orb216514)

| | |
|-------------------------|---|
| Description | Bradykinin Potentiator C peptide |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Protein Sequence | pGlu-Gly-Leu-Pro-Pro-Gly-Pro-Pro-Ile-Pro-Pro-OH |
| Formula | $C_{51}H_{77}N_{11}O_{13}$ |
| Purity | > 95% |
| Source | Synthetic |
| MW | 1052.26 |
| CAS Number | 30953-20-9 |
| Expiration Date | 12 months from date of receipt. |

