
Product Datasheet

Sheep IgG Rat IgG Fab (TRITC) antibody (orb21659)

Description

Tetramethylrhodamine isothiocyanate-

sheep antiserum to rat Fab...

Species/Host

Sheep

Conjugation

TRITC

Preservatives

TRITC-coupled purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added as it may interfere with he antibody activity.

Form/Appearance

TRITC-coupled purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added as it may interfere with he antibody activity.

Storage

Storage: The lyophilized conjugate is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the immunoconjugate

Note

For research use only

Clonality

Polyclonal

Source

Purified Fab of normal IgG isolated form rat serum. Freund's complete adjuvant is used in the first step of the immunization procedure.

Hazard Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.