

Product Datasheet

Gastrin Receptor I Human peptide (orb216604)

Description	Gastrin Receptor I Human peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Selective CCK2 agonist. Endogenous peptide, produced in the stomach and capable of stimulating gastric acid secretion.
Protein Sequence	pGlu-Gly-Pro-Trp-Leu-Glu-Glu-Glu-Glu-Glu-Ala-Tyr-Gly-Trp-Met-Asp-Phe-NH ₂
Formula	C ₉₇ H ₁₂₅ N ₂₀ O ₃₁ S
Purity	> 95%
Source	Synthetic
MW	2098.22
CAS Number	10047-33-3
Expiration Date	12 months from date of receipt.