


Product Datasheet

Hexarelin peptide (orb216643)

Description	Hexarelin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Highly potent growth hormone secretagogue receptor agonist (IC ₅₀ = 18 pM). Inhibits apoptosis. Ghrelin analogue. Cardioprotective agent.
Protein Sequence	His-D-2-Me-Trp-Ala-Trp-D-Phe-Lys-NH ₂
Formula	C ₄₇ H ₅₈ N ₁₂ O ₆
Purity	> 95%
Source	Synthetic
MW	887
CAS Number	140703-51-1
Expiration Date	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)