

Product Datasheet

Goserelin Acetate peptide (orb216725)

Description	Goserelin Acetate peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Goserelin peptide
Protein Sequence	Pyr-His-Trp-Ser-Tyr-D-Ser (tBu)-Leu-Arg-Pro-AzaGly-NH ₂
Formula	C ₅₉ H ₈₄ N ₁₈ O ₁₄
Purity	> 95%
Source	Synthetic
MW	1269.4
CAS Number	65807-02-5
Expiration Date	12 months from date of receipt.

