

Product Datasheet

ADARB1 peptide (orb216880)

Description	The peptide is used to block anti-ADAR2 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	ADARB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213524.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q91ZS8 , P78563 , P51400
Entrez	25367 , 104 , 110532
Expiration Date	12 months from date of receipt.