

---

## Product Datasheet

### Rabbit IgG Monkey IgG Fab (HRP) antibody (orb21699)

**Description**

Horseradish peroxidase-conjugated IgG

to monkey Fab of IgG. ...

**Species/Host**

Rabbit

**Reactivity**

Monkey

**Conjugation**

HRP

**Tested**

DOT, ELISA, ICC, IHC-P, WB

**Applications**
**Storage**

Storage: The lyophilized conjugate is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the immunoconjugate

**Note**

For research use only

**Clonality**

Polyclonal

**Source**

Purified normal Fab of polyclonal rhesus monkey IgG. Freund's complete adjuvant is used in the first step of the immunization procedure.

**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.

**Expiration Date**

12 months from date of receipt.