

## Product Datasheet

### CDC34 peptide (orb217059)

<b>Description</b>	The peptide is used to block anti-CDC34 antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Primate, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	CDC34
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb213703.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P49427</a> , <a href="#">Q8CFI2</a>
<b>Entrez</b>	<a href="#">997</a> , <a href="#">216150</a>
<b>Expiration Date</b>	12 months from date of receipt.