

Product Datasheet

CEBPD peptide (orb217081)

Description	The peptide is used to block anti-C/EBP delta/epsilon antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	CEBPD; CEBPE
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213725.
Purity	> 85%
Source	Synthetic
Uniprot ID	P56261 , Q6PZD9 , Q15744 , Q00322 , P49716 , Q03484
Entrez	1052 , 1053 , 25410 , 110794 , 12609 , 25695
Expiration Date	12 months from date of receipt.