

## Product Datasheet

### CREB1 peptide (orb217130)

<b>Description</b>	CREB1 peptide
<b>Reactivity</b>	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
<b>Tested Applications</b>	Blocking
<b>Target</b>	CREB1
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb213774.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q01147</a> , <a href="#">P15337</a> , <a href="#">P16220</a>
<b>Entrez</b>	<a href="#">1385</a> , <a href="#">12912</a> , <a href="#">81646</a>
<b>Expiration Date</b>	12 months from date of receipt.