

Product Datasheet

CYP4A11 peptide (orb217186)

| | |
|----------------------------|--|
| Description | The peptide is used to block anti-Cytochrome P450 4A11/22 antibody reactivity. |
| Reactivity | Human, Mouse, Rat |
| Tested Applications | Blocking |
| Target | CYP4A11; CYP4A22 |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb213830. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | Q02928 , Q5TCH4 |
| Entrez | 1579 , 284541 |
| Expiration Date | 12 months from date of receipt. |