

Product Datasheet

CYP11A1 peptide (orb217188)

Description	The peptide is used to block anti-Cytochrome P450 11A1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CYP11A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213832.
Purity	> 85%
Source	Synthetic
Uniprot ID	P14137 , P05108 , Q9QZ82
Entrez	13070 , 1583 , 29680
Expiration Date	12 months from date of receipt.