

Product Datasheet

HOXA11 peptide (orb217403)

Description	The peptide is used to block anti-HOXA11/D11 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	HOXA11; HOXD11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214047.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31270 , P31277 , P31311 , P23813
Entrez	15431 , 3237 , 3207 , 15396
Expiration Date	12 months from date of receipt.