
Product Datasheet

Rabbit IgG Human Ig Lambda Light Chain (FITC) antibody (orb21755)

Description

Fluorescein isothiocyanate-

polyclonal rabbit antiserum to free and bound ...

Species/Host

Rabbit

Reactivity

Human

Conjugation

FITC

Preservatives

Fluorochrome coupled purified rabbit IgG lyophilized from a solution in phosphate buffered saline (pH 7.2). No preservative added, as it may interfere with the antibody activity.

Form/Appearance

Fluorochrome coupled purified rabbit IgG lyophilized from a solution in phosphate buffered saline (pH 7.2). No preservative added, as it may interfere with the antibody activity.

Storage

Storage: The lyophilized product is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

Note

For research use only

Clonality

Polyclonal

Source

A pool of purified Bence Jones lambda proteins isolated from human urine. Freund's complete adjuvant is used in the first step of the immunization procedure.