

Product Datasheet

NKX3-1 peptide (orb217676)

Description	The peptide is used to block anti-NKX3-1 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	NKX3-1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214320.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99801
Entrez	4824
Expiration Date	12 months from date of receipt.