

Product Datasheet

CNTN3 peptide (orb217714)

Description	The peptide is used to block anti-Contactin 3 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CNTN3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214358.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62682 , Q07409 , Q9P232
Entrez	18488 , 54279 , 5067
Expiration Date	12 months from date of receipt.