

## Product Datasheet

### SCN8A peptide (orb217914)

<b>Description</b>	The peptide is used to block anti-Nav1.6 antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	SCN8A
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb214558.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q9UQD0</a> , <a href="#">Q9WTU3</a> , <a href="#">O88420</a>
<b>Entrez</b>	<a href="#">20273</a> , <a href="#">29710</a> , <a href="#">6334</a>
<b>Expiration Date</b>	12 months from date of receipt.