

Product Datasheet

SSX1 peptide (orb217970)

Description	The peptide is used to block anti-SSX-pan antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	SSX1; SSX2; SSX3; SSX4; SSX5; SSX6; SSX7; SSX8; SSX9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214614.
Purity	> 85%
Source	Synthetic
Uniprot ID	O60225 , Q99909 , Q7RTT6 , Q16384 , Q7RTT5 , Q7RTT4 , O60224 , Q16385 , Q7RTT3
Entrez	6756 , 727837 , 548313 , 280658 , 6759 , 10214 , 6757 , 6758
Expiration Date	12 months from date of receipt.