

Product Datasheet

CNTN2 peptide (orb217998)

| | |
|----------------------------|--|
| Description | The peptide is used to block anti-Contactin 2 antibody reactivity. |
| Reactivity | Human, Mouse, Rat |
| Tested Applications | Blocking |
| Target | CNTN2 |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb214642. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | Q61330 , Q02246 , P22063 |
| Entrez | 21367 , 6900 , 25356 |
| Expiration Date | 12 months from date of receipt. |