

Product Datasheet

UBE2D2 peptide (orb218067)

Description	The peptide is used to block anti-UBE2D2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	UBE2D2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214711.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62837 , P62839 , P62838
Entrez	641452 , 7322 , 56550
Expiration Date	12 months from date of receipt.