

Product Datasheet

UBE2I peptide (orb218068)

Description	The peptide is used to block anti-UBC9 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	UBE2I
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214712.
Purity	> 85%
Source	Synthetic
Uniprot ID	P63279 , P63281 , P63280
Entrez	25573 , 22196 , 7329 , 102641751
Expiration Date	12 months from date of receipt.