

Product Datasheet

USP11 peptide (orb218129)

Description	The peptide is used to block anti-USP11 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	USP11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214773.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99K46 , P51784
Entrez	236733 , 8237
Expiration Date	12 months from date of receipt.