

Product Datasheet

ADAMTS2 peptide (orb218260)

Description	ADAMTS2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ADAMTS2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214904.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8C9W3 , O95450
Entrez	216725 , 9509
Expiration Date	12 months from date of receipt.