

Product Datasheet

RNF144A peptide (orb218287)

Description	The peptide is used to block anti-RNF144A antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RNF144A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214931.
Purity	> 85%
Source	Synthetic
Uniprot ID	P50876
Entrez	9781
Expiration Date	12 months from date of receipt.