

Product Datasheet

OXSR1 peptide (orb218303)

Description	The peptide is used to block anti-OXSR1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	OXSR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214947.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6P9R2 , O95747
Entrez	108737 , 9943
Expiration Date	12 months from date of receipt.