

Product Datasheet

ABCA7 peptide (orb218334)

Description	The peptide is used to block anti-ABCA7 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ABCA7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214978.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8IZY2 , Q91V24 , Q7TNJ2
Entrez	10347 , 299609 , 27403
Expiration Date	12 months from date of receipt.