

Product Datasheet

USP19 peptide (orb218377)

Description	The peptide is used to block anti-USP19 antibody reactivity.
Reactivity	Human, Rabbit, Rat
Tested Applications	Blocking
Target	USP19
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215021.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6J1Y9 , O94966
Entrez	361190 , 10869
Expiration Date	12 months from date of receipt.