

Product Datasheet

41888 peptide (orb218430)

Description	The peptide is used to block anti-Septin 6 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SEPT6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215074.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9R1T4 , Q14141
Entrez	23157 , 56526
Expiration Date	12 months from date of receipt.