

## Product Datasheet

### IPPK peptide (orb218624)

<b>Description</b>	The peptide is used to block anti-IPPK antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Primate, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	IPPK
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb215268.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q6P1C1</a> , <a href="#">Q9H8X2</a> , <a href="#">Q5PXE9</a>
<b>Entrez</b>	<a href="#">306808</a> , <a href="#">64768</a> , <a href="#">75678</a>
<b>Expiration Date</b>	12 months from date of receipt.