

Product Datasheet

USP42 peptide (orb218675)

Description	The peptide is used to block anti-USP42 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	USP42
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215319.
Purity	> 85%
Source	Synthetic
Uniprot ID	B2RQC2 , D3ZU96 , Q9H9J4
Entrez	76800 , 288482 , 84132
Expiration Date	12 months from date of receipt.