

Product Datasheet

SARNP peptide (orb218676)

Description	The peptide is used to block anti-CIP29 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SARNP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215320.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9D1J3 , Q498U4 , P82979
Entrez	362819 , 66118 , 84324
Expiration Date	12 months from date of receipt.