

Product Datasheet

Mouse KLRG1/MAFA Antibody : PE [Out of stock] (orb2186816)

Description	Mouse KLRG1/MAFA Antibody : PE [Out of stock]
Species/Host	Hamster
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	IL-2-activated NK cells (A-LAK cells) from B6 (H-2b) mice
Storage	Store at 2-8C
Note	For research use only
Application notes	Application Info: Flow cytometry 1, Immunoprecipitation1. Dilution Info: Flow Cytometry: Purified antibodyFluorescein conjugateBiotin conjugateR-phycoerythrin conjugate6 cells6 cells6 cells6 cells
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q96E93 , O88713
NCBI	NP_058666.1
Expiration Date	12 months from date of receipt.