

Product Datasheet

CBLN4 peptide (orb218731)

Description	The peptide is used to block anti-Cerebellin 4 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	CBLN4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215375.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NTU7 , Q8BME9
Entrez	228942 , 140689
Expiration Date	12 months from date of receipt.