

Product Datasheet

FOXR1 peptide (orb218763)

Description	The peptide is used to block anti-FOXR1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FOXR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215407.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6PIV2 , Q3UTB7
Entrez	283150 , 382074
Expiration Date	12 months from date of receipt.