

Product Datasheet

PLSCR1 peptide (orb218860)

Description	The peptide is used to block anti-PLSCR1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	PLSCR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215504.
Purity	> 85%
Source	Synthetic
Uniprot ID	P58195 , Q9JJ00 , O15162
Entrez	117540 , 5359 , 22038
Expiration Date	12 months from date of receipt.