

Product Datasheet

PPP1R1A peptide (orb218863)

Description	The peptide is used to block anti-IPP1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PPP1R1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215507.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13522 , Q9ERT9 , P19103
Entrez	58977 , 58200 , 5502
Expiration Date	12 months from date of receipt.