

Product Datasheet

RAB14 peptide (orb218925)

Description	The peptide is used to block anti-RAB14 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	RAB14
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215569.
Purity	> 85%
Source	Synthetic
Uniprot ID	P61107 , P61106 , Q91V41
Entrez	68365 , 51552 , 94197
Expiration Date	12 months from date of receipt.