

## Product Datasheet

### MED9 peptide (orb218929)

<b>Description</b>	The peptide is used to block anti-MED9 antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Primate, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	MED9
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb215573.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q9NWA0</a> , <a href="#">Q8VCS6</a>
<b>Entrez</b>	<a href="#">55090</a> , <a href="#">192191</a>
<b>Expiration Date</b>	12 months from date of receipt.