

## Product Datasheet

# Sheep IgG F(c) Antibody - Texas Red Conjugated (orb2190210)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Sheep IgG F(c) Antibody - Texas Red Conjugated   |
| <b>Conjugation</b>         | Texas Red  |
| <b>Tested Applications</b> | FC, FLISA, IF  |
| <b>Immunogen</b>           | Sheep IgG F(c) fragment  |
| <b>Preservatives</b>       | Lyophilized  |
| <b>Form/Appearance</b>     | Lyophilized  |
| <b>Storage</b>             | 2°C to 8°C -20°C or -80°C  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500 |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 625  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |