

## Product Datasheet

# Mouse IgG/IgA/IgM (H&L) Antibody - Biotin Conjugated (orb2190644)

<b>Description</b>	Mouse IgG/IgA/IgM (H&L) Antibody - Biotin Conjugated
<b>Species/Host</b>	Goat
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Mouse IgG IgA and IgM whole molecule
<b>Preservatives</b>	Lyophilized
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	2°C to 8°C -20°C or -80°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: Immunoblotting, ELISA, immunohistochemistry, immunomicroscopy as well as other antibody based assays using streptavidin or avidin conjugates requiring lot-to-lot consistency.
<b>Clonality</b>	Polyclonal
<b>MW</b>	455
<b>Expiration Date</b>	12 months from date of receipt.